

Quality Control Log – Glucometer

1. Keep glucometer(s) operator’s manual with this log.
2. Per operator’s manual, perform control solution testing when:
 - Using the meter for the first time.
 - Opening a new bottle of test strips.
 - The test strip bottle cap has been left open for a while.
 - The meter has been dropped.
 - You suspect the meter or test strips are not working properly.
 - Results are repeatedly unexpected/does not match how patient feels.
 - Other: _____
 - Other: _____

Date	Time	Staff Initials	Test Strip Lot #	Low Control			High Control			Out of Range (✓)	Action Taken/Comments
				Lot #	Range	Result	Lot #	Range	Result		

Print Name: _____ Signature: _____ Initials: _____
 Print Name: _____ Signature: _____ Initials: _____
 Print Name: _____ Signature: _____ Initials: _____
 Print Name: _____ Signature: _____ Initials: _____