Sterilized Packages Inspection Log

Sterilized items are considered sterile until use unless an event causes contamination. Sterilized items are not considered sterile if the package is opened, wet or moist, discolored, or damaged. Compromised packages shall be removed from the sterile package storage area and immediately repackaged, relabeled, and re-sterilized. Entries below indicate routine evaluation of the integrity and condition of sterilized packages according to the site's written schedule. These inspection logs shall be accessible on site for at least three years.

Directions:

Year: _____

Staff shall initial the appropriate box (month and day) indicating when routine inspections occurred. Staff shall initial and sign the bottom of this form for proper identification.

	1	E.I.	1 14		1 14	1	11		6	0.1	.	.
	Jan	Feb	Mar	Apr	May	Jun	Jul	Aug	Sep	Oct	Nov	Dec
1												
2												
3												
3 4 5												
6												
7												
8												
9												
10												
11												
12												
13												
14												
15												
16												
17												
18												
19												
20 21												
21												
22												
23												
24												
25												
26												
27												
28												
29												
30												
31												
		•										

Print Staff Name: ______ Staff Signature: ______ Initials: _____

Print Staff Name: ______ Staff Signature: _____ Initials: ____

Print Staff Name: ______ Staff Signature: _____ Initials: ____

Print Staff Name: ______ Staff Signature: _____ Initials: ____

CABC-CD-050390-24	March 2024
C. IDC CD 000070 2 1	

Staff Signature / Initials: